

Assay of Clopidogrel by Using HPLTC Method

Dr. D. S. Ghotekar¹, Vishal N. Kushare²

¹Department of Chemistry, N.V.P.Mandal's, Arts, Commerce and Science College Lasalgaon

Tal-Niphad, Dist-Nashik, Maharashtra, India

²Department of Pharmaceutics, Maratha Vidya Prasarak Samaja's Institute of Pharmaceutical

Sciences, Adgaon, Nashik, Maharashtra, India

ABSTRACT

The present paper describes stability indicating high-performance thin-layer chromatography (HPTLC) assay method for clopidogrel in bulk drugs. The method employed TLC aluminium plates precoated with silica gel 60F-254 as the stationary phase. The solvent system consisted of toluene: methanol: triethylamine (6.5: 4.0: 0.1 v/v/v). The system was found to give compact spot for clopidogrel (R_f value of 0.40 ± 0.010). Densitometric analysis of Clopidogrel was carried out in the absorbance mode at 254 nm. The linear regression analysis data for the calibration plots showed good linear relationship with $r^2 = 0.888$ with respect to peak area in the concentration range 30 - 120 ng/spot.

Keywords : Clopidogrel, Validation, HPTLC

I. INTRODUCTION

High Performance Thin Layer Chromatography (HPTLC) is the most powerful advanced form of Thin Layer Chromatography (TLC) and consists of chromatographic layers of utmost separation efficiency and the application of sophisticated instrumentation for all steps in the procedure include accurate sample application, standardized reproducible chromatogram development and software controlled evaluation. HPTLC is a concept that includes a widely standardized methodology based on scientific facts as well as the use of validated methods for qualitative and quantitative analysis. HPTLC meets all quality requirements for today's analytical labs, to increase the resolution and to allow more accurate quantitative measurements.

HPTLC is the most advanced form of modern TLC. It uses HPTLC plates featuring small particles with a narrow size distribution which results in homogenous layers with a smooth surface to be obtained. HPTLC uses smaller plates (10×10 or 10×20 cm). HPTLC plates provide improved resolution, higher detection sensitivity, and improved in situ quantification and are used for industrial pharmaceutical densitometric quantitative analysis. Normal phase adsorption TLC on silica gel with a less polar mobile phase, such as chloroform– methanol, has been used for more than 90% of reported analysis of pharmaceuticals and drugs.

The selection of mobile phase is based on adsorbent material used as stationary phase and physical and chemical properties of analyte. The mobile-phase systems are used based on their diverse selectivity

properties are diethyl ether, methylene chloride, and chloroform combined individually or together with hexane as the strength adjusting solvent for normal-phase TLC and methanol, acetonitrile, and tetrahydrofuran mixed with water for strength adjustment in reversed-phase TLC. Separations by ion pairing on C-18 layers are done with a mobile phase such as methanol-0.1 M acetate buffer (pH 3.5) containing 25 mM sodium pentanesulfonate (15.5:4.5).

I. EXPERIMENTAL WORK

❖ Standard Solutions

• Solution A (Stock Standard Solution)

Accurately weighed quantity of Clopidogrel (10.0 mg) was dissolved in methanol to make 10.0 ml solution. (conc.: 1.0 mg/ml).

• Solution B (Working Standard Solution)

Accurately measured 1.0 ml of solution A was diluted to 100.0 ml with methanol (conc.: 10.0 µg/ml).

❖ Optimization of chromatographic conditions

• Optimization of mobile phase


Aliquot portions of working standard solution (5 µl) were applied on TLC plates. Various pure solvents with varying polarity and their mixtures were tried for optimum movement of drug with sharp symmetrical peak. After trying several permutations and combinations, the mobile phase containing toluene: methanol: triethylamine (6.5: 4.0: 0.1 v/v/v) was found to be most satisfactory as it gave sharp symmetrical peaks for the drug with R_F values 0.40 ± 0.010 .

• Optimization of wavelength for densitometric evaluation


Aliquot portion of working standard solution B (5 µl) was applied on TLC plate. The plate was developed using optimized mobile phase. After development, the plate was removed and dried with the help of hair dryer. The migrated band was scanned over the wavelength range 200- 400

nm in an absorbance/reflectance mode and an *in situ* UV-absorption spectrum of drug was obtained.

A 254.0 nm was selected as scanning wavelength as it gave maximum absorption for the drug.



(a)



(b)

Figure 1 : (a) HPTLC chromatogram and (b) *In situ* UV Spectrum of clopidogrel



• Chamber saturation time

The chamber saturation time was optimized by allowing the spotted plates to equilibrate for varying time with vapours of mobile phase in twin trough chamber and then developing them using optimized mobile phase. The optimum saturation time was found to be 15 min, which resulted in to dense and compact spot.



On the basis of exhaustive experimentation the chromatographic conditions optimised are as follows:

Stationary Phase	Pre-coated Silica Gel 60 F254 TLC Plate (10 x10 cm)
Thickness	200 μ m
Mobile Phase	Toluene Methanol
Mode of Application	Band
Band Width	5mm
Sample Volume	5 μ l
Separation Technique	Ascending
Development Chamber	Twin trough glass chamber , 10 x 10 cm
Chamber Saturation Time	15 min
Migration Distance	90 mm
Detection	UV Densitometric Scanning
Scanning Mode	Absorbance/ Reflectance
Scanning Wavelength	254.0 nm
Scanning Speed	20 mm/s
Slit Dimension	4.0 x 0.45 mm
Temperature	25 \pm 3 $^{\circ}$ C


The samples of stress studies were spotted (5 μ l bands) on TLC plates and developed and scanned under optimized conditions. The study has shown adequate resolution of parent drug and its degradation products from one another under optimized chromatographic conditions by normal phase mode. Hence these chromatographic conditions were finalized for further experimentation.




A) Acid (1M HCL, 24 h reflux) and Spectra of degradant




B) Base (1M NaOH, 24 h reflux) and Spectra of degradant



C) Neutral stress, 24h reflux



E) Photolytic Stress, (15 days, 4 lac lux h) F) Thermal Stress (dry heat, 15 days at 700c



D) Oxidative stress (3% H₂O₂, 48 h at R.T.)

Figure 2 : HPTLC chromatograms of forced degraded samples (A-F).






Figure 3 : Chromatogram of mixed degradation sample

❖ Study of Linearity of response

Aliquots portions of standard solution B (3-12 µl) were applied on TLC plate and chromatograms were developed under optimized chromatographic conditions. The linear regression curves are depicted along with correlation coefficient, slope and y-intercept. The curves were found to be linear between concentration ranges 30-120 ng/spot.



(a)



(b)

Figure 4 : Linearity (a) by area (b) by height

Concentration range	30- 120 ng/spot	
Parameter	Height	Area
Slope	1.176	79.862
Y-intercept	92.498	479.765
Correlation coefficient	0.888	0.888

Table 1 : Results of Linearity studies

REFERENCES

- [1]. Sethi, P. D., In HPTLC-Quantitative Analysis of Pharmaceutical Formulations, CBS Publisher and Distributor, New Delhi, IInd ed., 1996, p. 1-19.
- [2]. D.H. Shewiyoa,b,c, E. Kaaleb, P.G. Rishab, B. Dejaegherc, J. Smeyers-Verbekec, Y. Vander Heyden. HPTLC methods to assay active ingredients in pharmaceutical formulations: A

review of the method development and validation steps. Journal of Pharmaceutical and Biomedical Analysis, 2012, 66, p. 11– 23.

- [3]. Nadig, D. E., "Preparation of drug sample for Analysis", Handbook of Pharmaceutical Analysis; Ohannesian, L., New Jersey, 2002; p. 1-3.
- [4]. Sonia K, Beddi Bhavya shree, Dr.K.S.Lakshmi, "HPTLC Method Development and Validation: An Overview", J. Pharm. Sci. & Res. Vol. 9(5), 2017,652-657
- [5]. Vishal N. Kushare, Sachin S Kushare, Development of Validated Stability Indicating HPTLC Method for Assay of Ozagrel and its Pharmaceutical Formulations, International Journal of Scientific Research in Science, Engineering and Technology. November-December-2018 4 (11): 36-48]
- [6]. Vishal N. Kushare, Sachin S Kushare, Sagar V Ghotekar, Development of Validated UV Spectrophotometric Method for Assay ofOzagrel and its Pharmaceutical Formulations. International Journal of Scientific Research in Science, Engineering and Technology. November-December-2018 4 (11): 69-77]

Cite this article as :

Dr. D. S. Ghotekar, Vishal N. Kushare, " Assay of Clopidogrel by Using HPLTC Method", International Journal of Scientific Research in Science and Technology (IJSRST), Print ISSN : 2395-6011, Online ISSN : 2395-602X, Volume 4, Issue 11, pp.557-561, November-December-2018. Available at doi : <https://doi.org/10.32628/IJSRST12848252>
Journal URL : <http://ijsrst.com/IJSRST12848252>